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Material

- 1 [Acrylic sheet](#)

Step 1: Set your laser cutting machine settings

Laser-cutting acrylics is relatively simple. The basic steps for achieving good-quality cut acrylic parts are listed here.

Settings for laser cutting acrylics will vary according to the machine. The table below lists the settings we used:

Setting	Value
Material	Acrylic 3mm
Power (W)	40
Speed (mm/s)	20

We recommend testing the parameters for acrylic laser cutting using a test matrix. Download and laser cut the test matrix (page 44) file. This will only use about an 8x9x0.5cm acrylic sheet.

The result should look like this:

Step 2: Laser cutting

Now you have tested your laser cutting machine and acrylic sheet you can laser cut the following parts:

Remove film

- If there is protective film on the acrylic, remove it now.
- Your plates are now prepared and ready to use.

Use cases

Please Note: This site is under construction.

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
Troubleshooting

Please Note: This site is under construction.

Use your laser module

Please Note: This site is under construction.

Print the plastic parts

Tool

- 1 [Precision wire cutter](#)
- 1 [RepRap-style printer](#)
- 1 [Utility knife](#) - Not a scalpel!

Material

- 200 gr of [PLA filament](#)

Step 1: Set your printer settings

Almost all station parts can be printed out of PLA filament on most RepRap-style printers.

We recommend the following printer settings:

Setting	Value
Material	PLA
Material Colour	Black
Material Temperature	Recommended by the PLA brand
Layer height	0.2mm or less
Infill	Printer default or more
Brim	Recommended for all parts
Slice gap closing radius	0.001mm

Test whether your printer can print the pieces for this station or other open-source designs. Download and print the leg test file. This will only use about 5 grams of PLA.

The result should look like this (this has been printed with a brim):

As a general rule, strength is more important than surface finish, so very thin layers (less than 0.15mm or so) are unlikely to result in a station that performs any better, though it may improve the appearance.

Step 2: Printing

Now you have tested your 3D printer and filament you can print the following parts:

You can download all the STLs as a single zipfile.


Step 3: Clean-up of printed parts

Be careful when removing brim: To avoid injury, first remove the bulk of the brim without a knife. Remove the remaining brim with a peeling action as described below.

Carefully remove the printing brim from all parts. To remove brim:

1. Use precision wire cutters to remove most of the brim from the part.
2. Clean up the remaining brim with a utility knife:
 - Hold the knife in your dominant hand with 4 fingers curled around the handle, leaving your thumb free.
 - Hold the part in your other hand, as far away from the surface, to be cut as possible.
 - Support the part with the thumb of your dominant hand.
 - Place the blade on the surface to be cut, and carefully close your dominant hand, moving the blade, under control, towards your thumb.

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CERN Open Hardware Licence Version 2 - Weakly Reciprocal

Preamble

CERN has developed this licence to promote collaboration among hardware designers and to provide a legal tool which supports the freedom to use, study, modify, share and distribute hardware designs and products based on those designs. Version 2 of the CERN Open Hardware Licence comes in three variants: CERN-OHL-P (permissive); and two reciprocal licences: this licence, CERN-OHL-W (weakly reciprocal) and CERN-OHL-S (strongly reciprocal).

The CERN-OHL-W is copyright CERN 2020. Anyone is welcome to use it, in unmodified form only.

Use of this Licence does not imply any endorsement by CERN of any Licensor or their designs nor does it imply any involvement by CERN in their development.

1 Definitions

1.1 'Licence' means this CERN-OHL-W.

1.2 'Compatible Licence' means

- a) any earlier version of the CERN Open Hardware licence, or
- b) any version of the CERN-OHL-S or the CERN-OHL-W, or
- c) any licence which permits You to treat the Source to which it applies as licensed under CERN-OHL-S or CERN-OHL-W provided that on Conveyance of any such Source, or any associated Product You treat the Source in question as being licensed under CERN-OHL-S or CERN-OHL-W as appropriate.

1.3 'Source' means information such as design materials or digital code which can be applied to Make or test a Product or to prepare a Product for use, Conveyance or sale, regardless of its medium or how it is expressed. It may include Notices.

1.4 'Covered Source' means Source that is explicitly made available under this Licence.

1.5 'Product' means any device, component, work or physical object, whether in finished or intermediate form, arising from the use, application or processing of Covered Source.

1.6 'Make' means to create or configure something, whether by manufacture, assembly, compiling, loading or applying Covered Source or another Product or otherwise.

1.7 'Available Component' means any part, sub-assembly, library or code which:

- a) **is** licensed to You **as** Complete Source under a Compatible Licence; **or**
- b) **is** available, at the time a Product **or** the Source containing it **is** first Conveyed, to You **and** any other prospective licensees
 - i) with sufficient rights **and** information (including any configuration **and** programming files **and** information about its characteristics **and** interfaces) to enable it either to be Made itself, **or** to be sourced **and** used to Make the Product; **or**
 - ii) **as** part of the normal distribution of a **tool** used to design **or** Make the Product.

1.8 'External Material' means anything (including Source) which:

- a) **is only combined with Covered Source in such a way that it interfaces with the Covered Source using a documented interface which is described in the Covered Source; and**
- b) **is not a derivative of or contains Covered Source, or, if it is, it is solely to the extent necessary to facilitate such interfacing.**

1.9 'Complete Source' means the set of all Source necessary to Make a Product, in the preferred form for making modifications, including necessary installation and interfacing information both for the Product, and for any included Available Components. If the format is proprietary, it must also be made available in a format (if the proprietary tool can create it) which is viewable with a tool available to potential licensees and licensed under a licence approved by the Free Software Foundation or the Open Source Initiative. Complete Source need not include the Source of any Available Component, provided that You include in the Complete Source sufficient information to enable a recipient to Make or source and use the Available Component to Make the Product.

1.10 'Source Location' means a location where a Licensor has placed Covered Source, and which that Licensor reasonably believes will remain easily accessible for at least three years for anyone to obtain a digital copy.

1.11 'Notice' means copyright, acknowledgement and trademark notices, Source Location references, modification notices (subsection 3.3(b)) and all notices that refer to this Licence and to the disclaimer of warranties that are included in the Covered Source.

1.12 'Licensee' or 'You' means any person exercising rights under this Licence.

1.13 'Licensor' means a natural or legal person who creates or modifies Covered Source. A person may be a Licensee and a Licensor at the same time.

1.14 'Convey' means to communicate to the public or distribute.

2 Applicability

2.1 This Licence governs the use, copying, modification, Conveying of Covered Source and Products, and the Making of Products. By exercising any right granted under this Licence, You irrevocably accept these terms and conditions.

2.2 This Licence is granted by the Licensor directly to You, and shall apply worldwide and without limitation in time.

2.3 You shall not attempt to restrict by contract or otherwise the rights granted under this Licence to other Licensees.

2.4 This Licence is not intended to restrict fair use, fair dealing, or any other similar right.

3 Copying, Modifying and Conveying Covered Source

3.1 You may copy and Convey verbatim copies of Covered Source, in any medium, provided You retain all Notices.

3.2 You may modify Covered Source, other than Notices, provided that You irrevocably undertake to make that modified Covered Source available from a Source Location should You Convey a Product in circumstances where the recipient does not otherwise receive a copy of the modified Covered Source. In each case subsection 3.3 shall apply.

You may only delete Notices **if** they are no longer applicable to the corresponding Covered Source as modified by You and You may add additional Notices applicable to Your modifications.

3.3 You may Convey modified Covered Source (with the effect that You shall also become a Licensor) provided that You:

- a) retain Notices as required in subsection 3.2;
- b) add a Notice to the modified Covered Source stating that You have modified it, with the date and brief description of how You have modified it;
- c) add a Source Location Notice **for** the modified Covered Source **if** You Convey in circumstances where the recipient does not otherwise receive a copy of the modified Covered Source; **and**
- d) license the modified Covered Source under the terms and conditions of this Licence (or, as set out in subsection

8.3, a later version, **if** permitted by the licence of the original Covered Source). Such modified Covered Source must be licensed as a whole, but excluding Available Components contained in it or External Material to which it is interfaced, which remain licensed under their own applicable licences.

4 Making and Conveying Products

4.1 You may Make Products, and/or Convey them, provided that You either provide each recipient with a copy of the Complete Source or ensure that each recipient is notified of the Source Location of the Complete Source. That Complete Source includes Covered Source and You must accordingly satisfy Your obligations set out in subsection 3.3. If specified in a Notice, the Product must visibly and securely display the Source Location on it or its packaging or documentation in the manner specified in that Notice.

4.2 Where You Convey a Product which incorporates External Material, the Complete Source for that Product which You are required to provide under subsection 4.1 need not include any Source for the External Material.

4.3 You may license Products under terms of Your choice, provided that such terms do not restrict or attempt to restrict any recipients' rights under this Licence to the Covered Source.

5 Research and Development

You may Convey Covered Source, modified Covered Source or Products to a legal entity carrying out development, testing or quality assurance work on Your behalf provided that the work is performed on terms which prevent the entity from both using the Source or Products for its own internal purposes and Conveying the Source or Products or any modifications to them to any person other than You. Any modifications made by the entity shall be deemed to be made by You pursuant to subsection 3.2.

6 DISCLAIMER AND LIABILITY

6.1 **DISCLAIMER OF WARRANTY** -- The Covered Source and any Products are provided 'as is' and any express or implied warranties, including, but not limited to, implied warranties of merchantability, of satisfactory quality, non-infringement of third party rights, and fitness for a particular purpose or use are disclaimed in respect of any Source or Product to the maximum extent permitted by law. The Licensor makes no representation that any Source or Product does not or will not infringe any patent, copyright, trade secret or other proprietary right. The entire risk as to the use, quality, and performance of any Source or Product shall be with You and not the Licensor. This disclaimer of warranty is an essential part of this Licence and a condition for the grant of any rights granted under this Licence.

6.2 **EXCLUSION AND LIMITATION OF LIABILITY** -- The Licensor shall, to the maximum extent permitted by law, have no liability for direct, indirect, special, incidental, consequential, exemplary, punitive or other damages of

any character including, without limitation, procurement of substitute goods or services, loss of use, data or profits, or business interruption, however caused and on any theory of contract, warranty, tort (including negligence), product liability or otherwise, arising in any way in relation to the Covered Source, modified Covered Source and/or the Making or Conveyance of a Product, even if advised of the possibility of such damages, and You shall hold the Licensor(s) free and harmless from any liability, costs, damages, fees and expenses, including claims by third parties, in relation to such use.

7 Patents

7.1 Subject to the terms and conditions of this Licence, each Licensor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in subsections 7.2 and 8.4) patent licence to Make, have Made, use, offer to sell, sell, import, and otherwise transfer the Covered Source and Products, where such licence applies only to those patent claims licensable by such Licensor that are necessarily infringed by exercising rights under the Covered Source as Conveyed by that Licensor.

7.2 If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Covered Source or a Product constitutes direct or contributory patent infringement, or You seek any declaration that a patent licensed to You under this Licence is invalid or unenforceable then any rights granted to You under this Licence shall terminate as of the date such process is initiated.

8 General

8.1 If any provisions of this Licence are or subsequently become invalid or unenforceable for any reason, the remaining provisions shall remain effective.

8.2 You shall not use any of the name (including acronyms and abbreviations), image, or logo by which the Licensor or CERN is known, except where needed to comply with section 3, or where the use is otherwise allowed by law. Any such permitted use shall be factual and shall not be made so as to suggest any kind of endorsement or implication of involvement by the Licensor or its personnel.

8.3 CERN may publish updated versions and variants of this Licence which it considers to be in the spirit of this version, but may differ in detail to address new problems or concerns. New versions will be published with a unique version number and a variant identifier specifying the variant. If the Licensor has specified that a given variant applies to the Covered Source without specifying a version, You may treat that Covered Source as being released under any version of the CERN-OHL with that variant. If no variant is specified, the Covered Source shall be treated as being released under CERN-OHL-S. The Licensor may also specify that the Covered Source is subject to a specific version of the CERN-OHL or any later version in which case You may apply this or any later version of CERN-OHL with the same variant identifier published by CERN.

You may treat Covered Source licensed under CERN-OHL-W as licensed under CERN-OHL-S **if** and only **if** all Available Components referenced in the Covered Source comply with the corresponding definition of Available Component **for** CERN-OHL-S.

8.4 This Licence shall terminate with immediate effect if You fail to comply with any of its terms and conditions.

8.5 However, if You cease all breaches of this Licence, then Your Licence from any Licensor is reinstated unless such Licensor has terminated this Licence by giving You, while You remain in breach, a notice specifying the breach and requiring You to cure it within 30 days, and You have failed to come into compliance in all material respects by the end of the 30 day period. Should You repeat the breach after receipt of a cure notice and subsequent reinstatement, this Licence will terminate immediately and permanently. Section 6 shall continue to apply after any termination.

8.6 This Licence shall not be enforceable except by a Licensor acting as such, and third party beneficiary rights are specifically excluded. An error occurredPlease try again laterContact SupportClose Full-text Access

3D Models

List of Contents

1. Driver Mounting Plate
2. **A:** Used for two drivers model [INSERT MODEL]
3. **B:** Used for one driver model [INSERT MODEL A] and one driver model [INSERT MODEL B]
4. Power Supply Mounting Plate
5. Single Fiber Tray
6. Double Fiber Tray

Precision wire cutters

Precision wire cutters, sometimes called flush cutters, are wire cutters that come to a sharp point and have a flat back surface. This allows small wires (or 3D-printed supports) to be cut flush against another surface. These are very useful for removing small 3D-printed ties used to hold the microscope together during printing.

Utility knife

Any standard utility knife should work for simple jobs like removing a brim.

Any retractable or fixed blade utility knife should work. This is probably safer than scalpel, but should still be used with caution! A Stanley Titan fixed blade knife is a high quality option with replaceable blades that are fixed rather than retractable.

Rep-Rap style 3D printer

The stage has been designed to be printed on low-cost entry-level fused filament fabrication 3D printers. Many of these printers are Open-Hardware and originate from the Rep-Rap project.

A popular example is the Prusa i3 MK3S.

Laser cutting machine

High-precision tool that uses a focused laser beam to cut through or engrave various materials.

It offers versatility in cutting intricate shapes, patterns, and designs in materials like wood, acrylic, metal and fabric.


Laser cutting machines are used in industries like manufacturing, crafting, signage, and engineering for precise and efficient material processing.

PLA filament

Supplier: Esun

Test Matrix

Download test-matrix.pdf

Acrylic sheet

5mm thick plates.

Supplier: Novoacril

License

CERN Open Hardware Licence Version 2 - Weakly Reciprocal

Preamble

CERN has developed this licence to promote collaboration among hardware designers and to provide a legal tool which supports the freedom to use, study, modify, share and distribute hardware designs and products based on those designs. Version 2 of the CERN Open Hardware Licence comes in three variants: CERN-OHL-P (permissive); and two reciprocal licences: this licence, CERN-OHL-W (weakly reciprocal) and CERN-OHL-S (strongly reciprocal).

The CERN-OHL-W is copyright CERN 2020. Anyone is welcome to use it, in unmodified form only.

Use of this Licence does not imply any endorsement by CERN of any Licensor or their designs nor does it imply any involvement by CERN in their development.

1 Definitions

1.1 'Licence' means this CERN-OHL-W.

1.2 'Compatible Licence' means

a) any earlier version of the CERN Open Hardware licence, or

b) any version of the CERN-OHL-S or the CERN-OHL-W, or

c) any licence which permits You to treat the Source to which it applies as if it were licensed under CERN-OHL-S or CERN-OHL-W provided that on Conveyance of any such Source, or any associated Product You treat the Source in question as being licensed under CERN-OHL-S or CERN-OHL-W as appropriate.

1.3 'Source' means information such as design materials or digital code which can be applied to Make or test a Product or to prepare a Product for use, Conveyance or sale, regardless of its medium or how it is expressed. It must include Notices.

1.4 'Covered Source' means Source that is explicitly made available under this Licence.

1.5 'Product' means any device, component, work or physical object, whether in finished or intermediate form, arising from the use, application or production of Covered Source.

1.6 'Make' means to create or configure something, whether by manufacture

assembly, compiling, loading or applying Covered Source or another Product or otherwise.

1.7 'Available Component' means any part, sub-assembly, library or code w

a) is licensed to You as Complete Source under a Compatible Licence; or

b) is available, at the time a Product or the Source containing it is first Conveyed, to You and any other prospective licensees

i) with sufficient rights and information (including any configuration and programming files and information about its characteristics and interface to enable it either to be Made itself, or to be sourced and used to Make a Product; or

ii) as part of the normal distribution of a tool used to design or Make a Product.

1.8 'External Material' means anything (including Source) which:

a) is only combined with Covered Source in such a way that it interfaces with the Covered Source using a documented interface which is described in the Covered Source; and

b) is not a derivative of or contains Covered Source, or, if it is, it is solely to the extent necessary to facilitate such interfacing.

1.9 'Complete Source' means the set of all Source necessary to Make a Product in the preferred form for making modifications, including necessary installation and interfacing information both for the Product, and for any included Available Components. If the format is proprietary, it must also be made available in a format (if the proprietary tool can create it) which is viewable with a tool available to potential licensees and licensed under a licence approved by the Free Software Foundation or the Open Source Initiative. Complete Source need not include the Source of any Available Component, provided that You include in the Complete Source sufficient information to enable a recipient to Make or source and use the Available Component to Make the Product.

1.10 'Source Location' means a location where a Licensor has placed Covered Source, and which that Licensor reasonably believes will remain easily accessible for at least three years for anyone to obtain a digital copy.

1.11 'Notice' means copyright, acknowledgement and trademark notices, Source Location references, modification notices (subsection 3.3(b)) and all notices that refer to this Licence and to the disclaimer of warranties that are included in the Covered Source.

1.12 'Licensee' or 'You' means any person exercising rights under this Licence

1.13 'Licensor' means a natural or legal person who creates or modifies Covered Source. A person may be a Licensee and a Licensor at the same time.

1.14 'Convey' means to communicate to the public or distribute.

2 Applicability

2.1 This Licence governs the use, copying, modification, Conveying of Covered Source and Products, and the Making of Products. By exercising any right under this Licence, You irrevocably accept these terms and conditions.

2.2 This Licence is granted by the Licensor directly to You, and shall apply worldwide and without limitation in time.

2.3 You shall not attempt to restrict by contract or otherwise the rights granted under this Licence to other Licensees.

2.4 This Licence is not intended to restrict fair use, fair dealing, or any other similar right.

3 Copying, modifying and Conveying Covered Source

3.1 You may copy and Convey verbatim copies of Covered Source, in any medium provided You retain all Notices.

3.2 You may modify Covered Source, other than Notices, provided that You undertake to make that modified Covered Source available from a Source Location should You Convey a Product in circumstances where the recipient does not otherwise receive a copy of the modified Covered Source. In each case subsection 3.3 shall apply.

You may only delete Notices if they are no longer applicable to the corresponding Covered Source as modified by You and You may add additional Notices applicable to Your modifications.

3.3 You may Convey modified Covered Source (with the effect that You shall also become a Licensor) provided that You:

- a) retain Notices as required in subsection 3.2;
- b) add a Notice to the modified Covered Source stating that You have modified it, with the date and brief description of how You have modified it;
- c) add a Source Location Notice for the modified Covered Source if You Convey in circumstances where the recipient does not otherwise receive a copy of the modified Covered Source; and

d) license the modified Covered Source under the terms and conditions of Licence (or, as set out in subsection 8.3, a later version, if permitted the licence of the original Covered Source). Such modified Covered Source must be licensed as a whole, but excluding Available Components contained in it or External Material to which it is interfaced, which remain licensed under their own applicable licences.

4 Making and Conveying Products

4.1 You may Make Products, and/or Convey them, provided that You either provide each recipient with a copy of the Complete Source or ensure that each recipient is notified of the Source Location of the Complete Source. That Complete Source includes Covered Source and You must accordingly satisfy Your obligations set out in subsection 3.3. If specified in a Notice, the Product must visibly and securely display the Source Location on it or its packaging or documentation in the manner specified in that Notice.

4.2 Where You Convey a Product which incorporates External Material, the Source for that Product which You are required to provide under subsection 4.1 need not include any Source for the External Material.

4.3 You may license Products under terms of Your choice, provided that such terms do not restrict or attempt to restrict any recipients' rights under this Licence to the Covered Source.

5 Research and Development

You may Convey Covered Source, modified Covered Source or Products to a licensee entity carrying out development, testing or quality assurance work on Your behalf provided that the work is performed on terms which prevent the licensee from both using the Source or Products for its own internal purposes and disclosing the Source or Products or any modifications to them to any person other than You. Any modifications made by the entity shall be deemed to be made by You pursuant to subsection 3.2.

6 DISCLAIMER AND LIABILITY

6.1 **DISCLAIMER OF WARRANTY** -- The Covered Source and any Products are provided 'as is' and any express or implied warranties, including, but not limited to, implied warranties of merchantability, of satisfactory quality, non-infringement of third party rights, and fitness for a particular purpose or use are disclaimed in respect of any Source or Product to the maximum extent permitted by law. The Licensor makes no representation that any Source or Product does not or will not infringe any patent, copyright, trade secret or other proprietary right. The entire risk as to the use, quality, and performance of any Source or Product shall be with You and not the Licensor. This disclaimer of warranty is an essential part of this Licence and a condition for the grant of any rights granted under this Licence.

6.2 EXCLUSION AND LIMITATION OF LIABILITY -- The Licensor shall, to the extent permitted by law, have no liability for direct, indirect, special, incidental, consequential, exemplary, punitive or other damages of any kind, including, without limitation, procurement of substitute goods or services, loss of use, data or profits, or business interruption, however caused and on any theory of contract, warranty, tort (including negligence), product liability or otherwise, arising in any way in relation to the Covered Source, modified Covered Source and/or the Making or Conveyance of a Product, even if advised of the possibility of such damages, and You shall hold the Licensor free and harmless from any liability, costs, damages, fees and expenses, claims by third parties, in relation to such use.

7 Patents

7.1 Subject to the terms and conditions of this Licence, each Licensor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in subsections 7.2 and 8.4) patent license to Make, have Made, use, offer to sell, sell, import, and otherwise transfer the Covered Source and Products, where such licence applies only to those patent claims licensable by such Licensor that are necessarily infringed by exercising rights under the Covered Source as Conveyed by that Licensor.

7.2 If You institute patent litigation against any entity (including a counterclaim or counterclaim in a lawsuit) alleging that the Covered Source or a Product constitutes direct or contributory patent infringement, or You seek any other remedy that a patent licensed to You under this Licence is invalid or unenforceable, then any rights granted to You under this Licence shall terminate as of the date such process is initiated.

8 General

8.1 If any provisions of this Licence are or subsequently become invalid or unenforceable for any reason, the remaining provisions shall remain effective.

8.2 You shall not use any of the name (including acronyms and abbreviations), image, or logo by which the Licensor or CERN is known, except where needed to comply with section 3, or where the use is otherwise allowed by law. Any such permitted use shall be factual and shall not be made so as to suggest any kind of endorsement or implication of involvement by the Licensor or its personnel.

8.3 CERN may publish updated versions and variants of this Licence which it considers to be in the spirit of this version, but may differ in detail to address new problems or concerns. New versions will be published with a new version number and a variant identifier specifying the variant. If the Licence has specified that a given variant applies to the Covered Source without a version, You may treat that Covered Source as being released under any

of the CERN-OHL with that variant. If no variant is specified, the Covered Source shall be treated as being released under CERN-OHL-S. The Licensor also specify that the Covered Source is subject to a specific version of CERN-OHL or any later version in which case You may apply this or any later version of CERN-OHL with the same variant identifier published by CERN.

You may treat Covered Source licensed under CERN-OHL-W as licensed under CERN-OHL-S if and only if all Available Components referenced in the Covered Source comply with the corresponding definition of Available Component for CERN-OHL-S.

8.4 This Licence shall terminate with immediate effect if You fail to comply with any of its terms and conditions.

8.5 However, if You cease all breaches of this Licence, then Your Licence from any Licensor is reinstated unless such Licensor has terminated this Licence by giving You, while You remain in breach, a notice specifying the breach and requiring You to cure it within 30 days, and You have failed to come into compliance in all material respects by the end of the 30 day period. Should You repeat the breach after receipt of a cure notice and subsequent reinstatement, this Licence will terminate immediately and permanently. Section 6 shall continue to apply after any termination.

8.6 This Licence shall not be enforceable except by a Licensor acting as a licensor and third party beneficiary rights are specifically excluded.

Documentation powered by GitBuilding

©2023 Francisco Martinez, Matías Hurtado-Labarca, Pierre Padilla-Huamantínco, and Vicente Parot

Contact: librehub@uc.cl

Modular fiber-coupled laser illumination is released under CERN-OHL-W-2.0.